

## Appendix G | Label Symbols

Symbols may be found on the sensor, transmitter, and receiver package labels. These symbols tell you about the proper and safe use of the Dexcom G6. For a listing of what they mean, see the table below. You may also reference the Symbols Glossary at [dexcom.com/symbols](https://dexcom.com/symbols). This table shows what each symbol means.



Alternating Current



Batch/Lot Number



*Bluetooth*



Catalog Number



Caution




Class II Equipment



Consult Instructions for Use



Date of Manufacture



Direct Current



Do Not Reuse




Do Not Use if Package Is Damaged



Electrical Equipment Designed Primarily  
for Indoor Use



Humidity Limitation



Input

**IP22**

IP22: Protection Against Insertion of  
Large Objects and Dripping Water

**IP28**

IP28: Protection Against Insertion of  
Large Objects and Immersion in Water



Keep Away from Heat



Keep Dry



Manufacturer



MR Unsafe



Non-Ionizing Radiation

**PN**

Part Number

**Rx Only**

Prescription Required



Refer to Instruction Manual/Booklet

**SN**

Serial Number

**SB**

Ship-By Date

**STERILE R**

Sterile by Radiation



Temperature Limitation



Type BF Applied Part



Use By Date